



moshiri orthodontics
for a smile that lasts



AcceleDent®
aura



The AcceleDent® Aura with SoftPulse Technology™ is a device that has FDA approved to aid in tooth movement and may decrease the duration of your orthodontic treatment and can even help to make your orthodontics more comfortable. The Activator has an attached mouthpiece that is fitted around your existing appliances. With just 20 minutes per day, the AcceleDent® Aura generates small vibrations, or micropulses, enhancing the movements, allowing for faster tooth movement.

The AcceleDent® Aura with SoftPulse Technology™ has been demonstrated safe and reliable in U.S. clinical trial and has been FDA approved for use.

Appliance Usages:

- It is recommended that you use your **AcceleDent® Aura** for a continuous 20 minute period each day, as prescribed by Drs. Moshiri. After 10 minutes of initial wear, the activator will pulse three times. It will turn off after the 20 minutes has been completed.
- Bite pressure on the mouthpiece should be only firm enough to hold the appliance in place without using your hands. Avoid biting too forcibly, this will decrease the vibration forces of the teeth. The appliance does not need to be held once it's in place inside the mouth, but fee free to adjust the external unit with your hand occasionally.
- The Activator should be recharged every three to four days, when the LED indicator is quickly flashing orange, by using the wall charger or USB. **DO NOT WEAR** your **AcceleDent® Aura**, while charging. Avoid letting the battery completely discharge as this may affect the data reporting of the FastTrac Usage Report.
- The Activator is **NOT** water-proof. Do not use the **AcceleDent® Aura** while showering.
- The mouthpiece should be **disconnected** from the Activator before cleaning. The mouthpiece can be cleaned by rinsing in lukewarm water after each use. You may occasionally use a soft liquid dish soap, mild soap or toothpaste to gently brush the mouthpiece. **REMEMBER THE ACTIVATOR IS NOT WATER-PROOF**, so remove the mouthpiece from the Activator before cleaning.
- Keep your **AcceleDent® Aura** away from pets, especially dogs. They can damage or break the appliance. The **AcceleDent® Aura** is FDA cleared as a single-user prescription device.
- Contraindications for use: If patient is using osteoporosis drugs, has poor oral hygiene or has active periodontal disease that is not under full control for at least 3-4 months prior to the start of treatment.

Be sure and register your device and follow the Directions for Use manual, or if you have any questions, please call us at **(314) 4-BRACES**